EN

Distribution

P&R Medical
High Tech Science Park 1/13
University Hasselt 3590 B-Diepenbeek
Belgium
Tel: 0032(0) 11 23 46 81
Fax: 0032(0)11 23 46 82

Manufacturer: BHL-Medical
Quinta da Mata, Sete Casas
2670-350 Loures, Portugal 

REF: PR-20199 Transflux™ CT 20
PR-21379 Transflux™ CT 45
PR-21393 Transflux™ CT 100
PR-21410 Transflux™ CT 150
PR-20201 Transflux™ CT 60
PR-20202 Transflux™ CT 120
PR-20203 Transflux™ CT 250
PR-20206 Transflux™ MRI 150

Transflux™ application guide

1. Purpose

Transflux™ is a sterile single used product to be applied for CT, MRI and other procedures in radiology.

Transflux™ is to be used in combination with contrast injectors which have an automatic lock for pressures up to maximum 300 psi (20.615 bar or 2068.5 kPa), a Transset™ system and contrast media.


For further information regarding the injector, Transset and contrast media please consult the specific instructions for use.

2. Description of the Transflux™ system


Transflux™ is a patient line with one-way valves, safety zones, Clip, 3-LL, Reflux control system. The one-way valves offer the security of an anti contamination barrier between the patient and the transfer system.

3. Application


These systems are designed only for single use (one-way) to be used with a power injector up to a maximum pressure of 300 psi (20.615 bar or 2068.5 kPa) in radiology.

1. Open the peel package, take out the Transflux™ and remove the protective caps.
2. Keep Transflux™ vertically and fill the saline emulsion into the Transflux™ either manually with the syringe or automatically with the injector in order to remove any air bubbles.

3. Insert the needle into the vein or into a port-a-cath and make sure that the Luer fittings are firmly connected.
4. Connect Transflux™ with the needle.

5. Avoid subcutaneous injection or extravasations. Briefly disconnect the connection at the one-way valve proximal to the patient to check the blood reflux.
6. Connect Transflux™ with the transfer set and the injector. Briefly disconnect the connection at the one-way valve proximal to the patient to check the blood reflux. Prime the connector tube and make certain that the air is expelled.

AFTER INJECTION. Disconnect the patient line and leave the proximal part of the Transflux™ connected to the Transset™ as a safety cap. Close the clamp. Remove this part before installing the new Transflux™ for a new patient.

4. Contraindications

Not elsewhere specified.

5. Warnings symbols

Transflux™ is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Not for reuse
Any reuse of this product may cause a change in its mechanical and biological characteristics and could lead to product failure, allergic reactions or bacterial infections.

Do not use when package is damaged!
Patient or operator injury may result if package is opened or damaged or if damaged components are used. For such reason, it is very important that you visually inspect container and package before each use.

Do not re-sterilize the product!

- Air embolization can cause death or serious injury to the patient. Do not connect the patient to the injector until all trapped air has been cleared from the syringe and fluid path. Carefully read all instructions in the application guide to reduce the chance of air embolism.

- Do never exceed the indicated pressure of 300 psi (20.615 bar or 2068.5 kPa). Use of greater pressure or occlusions in the fluid path might result into leaks or ruptures of the Transflux.

- Patient or operator injury may result from contrast media leaks or tubing ruptures. Please ensure that the fluid path is open.

6. Precautions

Please take care that all connections are properly installed and assure they are all firm and secure. Please do not overtighten! This will help minimize leaks, disconnections and system damage.

7. Adverse effects

Not elsewhere specified


8. Storage

The product should be stored in its original package at a cool and dry place under normal temperature. Avoid direct sunlight.

9. Disclaimer

The manufacturer will not be liable for any damage caused by reuse of the product, abuse, improper handling and non-compliance with warning and handling instructions.

Identification of the graphic symbols

-  Caution! Consult accompanying documents
-  Consult instructions for use!
-  Use by YYYY-MM
-  Batch code
-  Order number
-  Do not use when package is damaged
-  Do not reuse!
-  Sterilized using ethylene oxide
-  Do not re-sterilize the product
-  No DEHP
-  Keep away from sunlight
-  Keep dry

Devices are worldwide patented and have Trade Mark.
*In certain jurisdictions based upon PCT/EP2006/02554, PCT/EP2009/06363 and PCT/EP2011/062179

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