	Vertrieb: P&R Medical High Tech Science Park 1/13 Universität Hasselt 3590 Diepenbeek – Belgien Tel.: +32(0)11234681 Fax: +32(0)11234682		Hersteller: BHL-Medical Quinta Da Mata. Sete-Casas. 2670-350 Loures, Portugal Fon +351 219821694 Fax +351 219821695	
				

	Distributed by P&R Medical High Tech Science Park 1/13 University Hasselt 3590 Diepenbeek – Belgium Tel: +32(0)11234681 Fax: +32(0)11234682		Manufacturer: BHL-Medical Quinta Da Mata. Sete-Casas. 2670-350 Loures, Portugal Fon +351 219821694 Fax +351 219821695	
				

	Distribuido por: P&R Medical High TECH Science Park 1/13 University Hasselt 3590 Diepenbeek – Bélgica Tel.: +32(0)11234681 Fax: +32(0)11234682		Fabricante: BHL-Medical Quinta Da Mata. Sete-Casas. 2670-350 Loures, Portugal Fon +351 219821694 Fax +351 219821695	
				

	Distribuido por: P&R Medical High TECH Science Park 1/13 University Hasselt 3590 Diepenbeek – Bélgica Tel.: +32(0)11234681 Fax: +32(0)11234682		Fabricado por: BHL-Medical Quinta Da Mata. Sete-Casas. 2670-350 Loures, Portugal Fon +351 219821694 Fax +351 219821695	
				

	Distribuido por: P&R Medical High TECH Science Park 1/13 University Hasselt 3590 Diepenbeek – Bélgica Tel.: +32(0)11234681 Fax: +32(0)11234682		Fabricado por: BHL-Medical Quinta Da Mata. Sete-Casas. 2670-350 Loures, Portugal Fon +351 219821694 Fax +351 219821695	
				

	Distribuido por: P&R Medical High TECH Science Park 1/13 University Hasselt 3590 Diepenbeek – Bélgica Tel.: +32(0)11234681 Fax: +32(0)11234682		Fabricado por: BHL-Medical Quinta Da Mata. Sete-Casas. 2670-350 Loures, Portugal Fon +351 219821694 Fax +351 219821695	
				

	Distribuido por: P&R Medical High TECH Science Park 1/13 University Hasselt 3590 Diepenbeek – Bélgica Tel.: +32(0)11234681 Fax: +32(0)11234682		Fabricado por: BHL-Medical Quinta Da Mata. Sete-Casas. 2670-350 Loures, Portugal Fon +351 219821694 Fax +351 219821695	
				

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

Suitable for reuse over a period of twelve consecutive hours.

The Transset contrast delivery system has been developed for delivering contrast media from a reservoir to the Transflux™ patient-line which is equipped with one-way valves and safety zones. The purpose of the system is to avoid the need to change power syringes for each patient procedure by eliminating the risk of contamination of the syringes (study: "Investigative Radiology – Volume 47, Number 4 "Study on the Microbial Safety of an Infusion Set for Contrast-Enhanced Imaging"). The benefit is a considerable reduction in cost, time and waste. At one end the Transset filling-and injection set is attached to the syringes of a power injector and to reservoirs of saline and contrast media. At the other end it is attached to a Transflux™ patient-line with one-way valves and safety zones. The combination of one-way valves allows filling the syringes and automatically draining the contrast and saline to the Transflux™ patient-line. This Transflux™ patient-line protects the Transset Filling & Injection set against contamination and consequently the Transset can be reused if the guidelines and instructions for use as described below are followed. To avoid contamination between the disconnection of the patient-lines.

1. Sterility of the Transset is maintained by disconnecting the distal part of the Transflux™ patient-line from the valve at the safety zone. (The so called disconnect feature).
2. Or protect the distal end of the Transset the sterile cap provided.
3. Open the peel-open pouch.
4. Connect the drip chambers to contrast and saline reservoirs.
5. Aspirate both syringes of the injector.
6. Flush the Transset system with normal saline.
7. Attach the set to the Transflux patient-line.
8. Ensure there is a permanent flow of the saline.
9. Execute the contrast injection.
10. After injection, disconnect the distal part of the Transflux™, patient-line leaving the safety zone connected until the next procedure as a protective cap. Alternatively, protect the distal end of the Transset with the sterile cap provided.
11. Start the cycle again.

For further information or instructions for the injector, the Transflux™ and the contrast material, please consult the specific instructions for use and 3D-animation and video on www.transflux.be

	Consult instructions for use !		No DEHP
	Keep dry		Max. use 12 hours
	Keep away from sunlight		Use by YYYY-MM
	Do not resterilize the product		Sterilized using ethylene oxide
	Do not use when package is damaged		Batch code
	Caution ! Consult accompanying documents		Order number

Patented

"Devices are worldwide patented" and have Trade Market.
"In certain jurisdictions based upon PCT/EP2006/062554, PCT/EP2009/066383 and PCT/EP2011/062179

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